

Optimal choices

Kem-En-Tec Diagnostics A/S,
Kuldysen 10, DK-2630 Taastrup
Tel.: +45 3927 1771, fax: +45 3920 0178

E-mail: support@kem-en-tec.com
http://www.kem-en-tec.com

11/09

SUBSTRATES	ENZ.		METHOD			SPECIFICATIONS					
	Horseradish Peroxidase	Alkaline Phosphatase	ELISA	Blotting - w: Western	Immunohistochemistry	Liquid, Ready-to-use	Fizzing, Buffer-containing tablets	Dissolvable in H ₂ O	Standard Tablets. Dissolved in Buffer	Staining	Cat. Number
AEC	X			W	X		5 mg		Red		4190
BCIP/NBT		X		W		X			Purple		4410
DAB, pH 5.0	X			W			10 mg		Brown		4150
DAB, pH 7.0	X				X		10 mg		Brown		4170
Fast Red		X		W	X				Red		4210
OPD	X		X				2 mg	2 or 10 mg	Yellow		4090/4110/4260
pNPP & pNPP ONE		X	X			X			Yellow		4400/4401
TMB ONE (& pre-stained)	X		X			X			Blue		4380/4430
TMB PLUS/TMB PLUS2 (& pre-st)	X		X			X			Blue		4390/4440//4395/4445
TMB X-tra (& pre-stained)	X		X			X			Blue		4800/4810
TMB Blotting/Membrane	X			W		X			Blue		4450/4460/4600
Buffers											
Antibody Enhancer	X	X	X	X	X	X					4660
HRP-StabilPLUS/HRP-StabilPLUS2	X		X	X	X	X					4530/4730
AP-StabilPLUS		X	X	X	X	X					4540
Sample PLUS2	X	X	X	X		X					4781
Synthetic Blocking Buffers	X	X	X	X		X					4520/4650
WellChampion	X	X		X		X					4900
Membrane Buffers	X	X		X		X					4750/4760